Hand over Analysis Results

With the help of VDDS-media communication, it is possible to transfer analysis results from the evaluation modules (2D and 3D) to the practice management. VDDS-media is defined and maintained by [Verband Deutscher Dentalsoftwarehersteller, VDDS](#).

The prerequisite is that their practice management system (PVS) supports the VDDS info import according to the VDDS-media standard. See [VDDS-media: PVS level 6 and MMOINFIMPORT](#).

Then OnyxCeph^{3™} can be configured to allow the transfer of analysis results to the PVS.

In the OnyxCeph^{3™} system options you should select your practice management system. For this purpose, it must be available as PMS with a MMOINFIMPORT entry in the VDDS_MMI.ini.

In the evaluation modules (2D and 3D), there should then be a button with the name of your PMS, with which the transfer of the current analysis can be triggered.

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```
[ONYXCEPH3_DATA]
TYPE=FRS                                     //Image type
TYPEID=0                                       //Image type index (the image
type name may change due to translation)
NAME=Österr. Ges. f. KFO                      //Name of analyse
DESCR=FRS-Analyse d. Österr. Gesellschaft f. KFO //Analysis long name
CODE=52694899                                    //Analysis ID (the analysis name
may change due to translation)

DATE=20100425                                    //Date of acquisition YYYYMMDD

COUNT=23                                         //Number of analysis results

V1_Name=1. Skelettale Diagnose                //1. Variable (Headline,
therefore no linked values)

V2_Name=NSAr                                      //2. Variable: Name
V2_Value=124.12                                    //2. Variable: Value (incl. 1
decimal)
V2_Unit=°                                         //2. Variable: Unit
V2_Code=0N3747X90                                //2. Variable: Unique ID (the
variable name may change due to translation)
V2_Text=124,1°                                     //2. Variable: Text-formatted
value (decimal separator defined by Windows-OS)
V3_Name=SArGo                                     //3. Variable
V3_Value=142.50

...
V6_Name=WT (Sum)                                 //6. Variable: A verbas result
V6_Value=horizontale Wachstumstendenz
V6_Unit=
V6_Code=0N3752X6750
...
```

The output is always in the language selected on the workstation.

Then the VDDS info import server of the PMS is opened with the communication file as parameter.

After the PMS has read the data, "READY" should be set to 1.

OnyxCeph^{3™} will delete all communication files at the latest at the next restart.

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