


Info Menu



The info menu can be found top right in the main window and provides direct calls for the following info and support options:

- Notifications *
- License information *
- Release information *
- Support tools *
- Client language setting

Entries marked by * require internet access. Right click into the info menu panel background opens the Proxy server settings window.

Button [Notifications]

	Notifications (Count)
	Notifications (Count) - Critical hint

► The following notifications are displayed:

- new release available
- expiring/expired license
- low memory (separated for client and server)
- Windows™ version outdated (separated for client and server)

If the notification list contains a critical message (e.g., too little memory), the icon turns red and the list of notifications is initially expanded.

Button [License]

	License active
	License will expire in less than 2 weeks
	License expired

► Left click provides context menu to open program [OnyxRegister](#) and to display the license information in a browser window. The license status can be displayed with [CTRL]+[SHIFT]+left click on menu item [License information].

Button [Release]

	Release up-to-date
	Release outdated - update required!
	No online access to update server

► On left click, a hint comes up stating that the update has to be performed on the server. Left click + [SHIFT][CTRL][ALT] provides a download option for the update in ZIP format. This download option on the client makes sense for users who don't have internet access on the server to run a direct online update.

Button [Support]

	OnyxSupport (Teamviewer) can be launched
	OnyxSupport is running or cannot be launched since another teamviewer instance is already running

► Left click opens the OnyxSupport window with the Support ID to be submitted. The display of personal data can be prevented via the context menu [▼].

In case another teamviewer instance is already running, this has to be closed first.



Button [Language]

	Language selection German
	Language selection English
	Language selection Russian
	Language selection Spanish
	Language selection Italian
	Language selection French
	Language selection Czech
	Language selection Hungarian
	Language selection Polish
	Language selection Portugal

► Left click opens the selection list with all available client languages.

To apply the selected client language, OnyxCeph has to be re-opened on the workstation.

Note: To change the language for cross-system terms and names (e.g. analysis names) the system language needs to be changed. This can only be done by program OnyxDBServer on the server.

From:
<http://onyxwiki.net/> - [OnyxCeph^{3TM} Wiki]
Permanent link:
http://onyxwiki.net/doku.php?id=en:info_menu
Last update: 2023/08/15 08:54



