


## Panel Object List

For each 3D scene (e.g., in 3D modules or in 3D full screens) an object list is available. Therein, all 3D subobjects of the used 3D dataset are listed. Single or multiple list entries (subobjects) can be selected by a left mouse click resp by [CTRL] + left mouse click.

If no subobject is selected, all actions, manipulations or other changes performed in panel [object options](#) or in panel [view options](#) are applied to the whole scene (all subobjects).

If single or multiple subobject are selected, all actions, manipulations or other changes performed in panel object options or in panel view options are applied to the selected subobjs only.




### Context menu:

In the object list context menu, several items for preselection, visibility, mesh type and representation, colorization etc. are available. For different modules, the context menu items can slightly differ depending on the characteristics of the module itself.

Entry
01 Select All
02 Unselect All
03 Invert selection
04 Hide Selected
05 Show Selected Only
06 Show All
07 Reset Position
08 Merge Vertices
09 Split Vertices
10 Colorize Regions

<b>Entry</b>	
<b>11</b>	Hide Coloration
<b>12</b>	Split Object
<b>13</b>	Merge Objects
<b>14</b>	Web Export
<b>15</b>	Hide Maxilla
<b>16</b>	Hide Mandibula
<b>17</b>	Show Maxilla Only
<b>18</b>	Show Mandibula Only
<b>19</b>	Plane Maxilla
<b>20</b>	Plane Mandibula
<b>21</b>	Mark All Teeth
<b>22</b>	Only Borer Teeth Visible
<b>23</b>	Shift Mandible
<b>24</b>	Select Brackets
<b>25</b>	Select Kylix
<b>26</b>	Select Base
<b>27</b>	Show Crowns Only
<b>28</b>	Show malocclusion Only
<b>29</b>	Show All


In Module [Adjust Image 3D](#) (only there) on top of the object list an icon menu appears err scheint which allows several actions to modify the object list entries. The icon bar replaces the former available panel Object Modification.



<b>Icon</b>	
<b>01</b>	Delete List Entry/Object
<b>02</b>	Rename List Entry/Object
<b>03</b>	Move List Entry/Object up
<b>04</b>	Move List Entry/Object down
<b>05</b>	Add List Entry/Object

## Mesh Inspection

By doubleclicking the statusbar at the bottom end of the panel, several topological inspection results for the single subobject surface mesh selected in the object list will be displayed → #3268.



From:

<http://onyxwiki.net/> - [OnyxCeph<sup>3™</sup> Wiki]



Permanent link:

[http://onyxwiki.net/doku.php?id=en:panel\\_objectlist](http://onyxwiki.net/doku.php?id=en:panel_objectlist)

Last update: **2016/07/12 10:24**