

Module Image Comparison 2D



Program version	BASIC	2D PRO	3D PRO	LAB	OMS
Availability	✗	■	■	■	■


✗ not included
 ■ included by
 default
 ■ optionally
 included

In module Image Comparison 2D, two or more findings of the same image type can be compared numerically and graphically. The numerical comparison is provided as a tabular comparison of the cephalometric measurement values of the selected analysis with or without clinical norms. The graphic comparison can be done as a superimposition of the tracing, as a segmented overlay and as a superimposition of harmony boxes. The tracings can be superimposed at reference points or (by means of gray value correlation) on growth-stable skeletal structures.

First Steps

1. Open module with 2 or more findings of the same image type
2. Configure numerical comparison
3. Select type and parameters for graphical superimposition
4. Adjust graphical superimposition

Module user interface



Links



From:
[https://www.onyxwiki.net/ - \[OnyxCeph3™ Wiki\]](https://www.onyxwiki.net/)

Permanent link:
https://www.onyxwiki.net/doku.php?id=en:image_comparison_2d&rev=1644220009

Last update: **2022/02/07 08:46**

