

Interface Configuration Digital Xray: DBSwin

Digital X-Ray

EasyDent

Carestream

Scanora

Orthometric/3DentStudio

VisiQuick

Vatech

VixWin

VixWin Scan

Sidexis

Orex

Dimaxis


Digora

CliniView

DBSwin

VistaScan

Romexis



☐ Duerr Dental DBSwin - Support Activated

DBSwin - Program Path

C:\Temp\xray\dbswin.exe

☒ Show Preview Images

☐ Smoothing

DBSwin Database

Test

Patient ID:

☒

Search

Patient:

☐

Images:

☐

Transfer Onyx Master Data

Transfer File:

Help

OK

Cancel

Parameter	Setting
[Support Activated]	Activates import icon button in module Add Image 2D
Programmpfad [...]	Lokaler Pfad zur ausführbaren Datei ..\DBSwin.exe
[Vorschaubilder anzeigen]	Aktivierung Thumbnail-Vorschau im Vorschaufenster Modul Bild hinzufügen 2D
[Glättung]	Verbesserung der Anzeige hochaufgelöster Bilder durch Resampling
[Test]	Testing database access
Übergabedatei [...]	Stammdatenübergabe mit Aufruf über Kommunikationsdatei]

Links



https://www.image-instruments.de/Onyx3/Bildimport_OnyxCeph3_aus_DBSWin.pdf

From:

<https://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

<https://www.onyxwiki.net/doku.php?id=en:dbswin>

Last update: **2022/06/27 07:33**

