

Interface Settings CT|PACS|KIS: HL7

CT / PACS / KIS

PACS HL7

☐ **HL7 - Support Activated**

Facility

Character Set

☐ **Import**

☐ HL7 Data Folder

☐ MLLP Server

Port

☐ **Export**

Application

Facility

Mode ☒ HL7 Data Folder ☐ MLLP Client

Path

☒ Embed document in message

Document Path

Help OK Cancel

Parameter	Setting
HL7 - Support Activated	Enable HL7 support here
Facility	Name of own HL7 facility
Character Set	Name of the character encoding used in the communication
Section Import	
Checkbox [HL7 Datenpfad]	OnyxCeph accepts messages via a communication directory
Checkbox [MLLP-Server]	OnyxCeph accepts messages via a MLLP server
Port	Receiving port of the OnyxCeph MLLP server
Section Export	
Application	Name of target application
Facility	Name of the receiving HL7 facility

Parameter	Setting
Mode	Send data via an HL7 communication directory <u>or</u> to an MLLP server
Path	Communication directory
Host	Address of the destination MLLP server
Port	Port of the destination MLLP server
Checkbox [Enable document in message]	Documents are transmitted as part of the message
Document Path	Directory for the transfer of documents

From:

<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

<http://onyxwiki.net/doku.php?id=en:hl7>

Last update: **2022/08/24 12:07**

