

Launch the archive process by hidden short cuts CTRL+doubleclick on label "Image Folder" in the OnyxDBServer window:

Archive Patient

Mark patient record for next archiving run or undo archiving this way to recover the archived record. Patient record will stay in the patient list if display of archived patient records is activated

Settings Menu:

- ☒ Show Archived Patients
- ☐ Sort Active Column
- ☒ Narrow Down Search Results
- ☐ Live Request
- ☒ Allow Combined Markers
- ☐ Cache List of Patients
- ☒ Update Cache Automatically
- ☐ Memorize Search Items
- ☒ Fit Columns

Patient's ID	Last Name	First Name	Gender	Model
10000	Demo	Daniela	Female	MN 550
10101	Demo	Daniela	Female	
10228	Demo	Daniela	Female	Doku
10432	Demo	Daniela	Female	
10000A	Demo	Daniela 이후 변	Female	MN 550
10240	Demo	Leipzig	Male	

From:

<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

http://onyxwiki.net/doku.php?id=en:img_4243_3

Last update: 2021/09/16 17:27

