

OnyxStat

Mit Hilfe des Dienstprogramms OnyxStat können Onyx-Datenbanken statistisch nach unterschiedlichen Kriterien durchsucht und ausgewertet werden. Die Abfrageergebnisse können gespeichert und im CSV-Format exportiert werden.

OnyxStat erfordert eine separate, einmalig kostenpflichtige Feischaltung.

Der Programmstart erfolgt aus der Programmgruppe OnyxCeph³™ im Windows-Startmenü am lokalen Arbeitsplatz.

Bedienoberfläche

The screenshot displays the OnyxStat application window. The main area shows 'Query 1' with a table of patient data. A secondary window, 'Query 2', is open, showing a detailed view of a specific patient's data. The 'Query 2' window includes a table with columns for Patient ID, Client, Name, First Name, Gender, Birth, Finding Date, Image Type, 2D/3D, Attributes, and Status. The table lists various findings and observations for patient 70020. On the right side of the interface, there is a 'View' panel with options for 'Color', 'Categories as Text', 'Sort By', 'Patient Filter', 'Data Filter', and 'Columns'. The 'Columns' list includes fields like Patient ID, Patient, Client, Name, First Name, Gender, Birth, Attributes, Finding Date, Image Type, 2D/3D, Status, Age, and GHA. At the bottom right, there is a 'Go to' section with an 'Image' thumbnail showing a lateral cephalogram.

Patient ID	Client	Name	First Name	Gender	Birth	Finding Date	Image Type	2D/3D	Attributes	Status
70020	Image Instru...	1	2	Male	24.04.2011	05.11.2012	Lateral Ceph	2D		Initial
7654321	Image Instru...	1	Patient(en)in...	Male	20.03.2003					
10015	Image Instru...	15								
10017	Image Instru...	17								
54321	Image Instru...	2								
10091	Image Instru...	2								
2000	Image Instru...	2D Ceph A								
ID001344	Image Instru...	3.7.1117.								
70016	Image Instru...	3D								
3Q	Image Instru...	3Q Techno								
10128	Image Instru...	3d								
10178	Image Instru...	3shape								
10118	Image Instru...	<SoX								
10115	Image Instru...	A								
AAAAA	Image Instru...	AA								
AAAAA	Image Instru...	AA								
AAAAA	Image Instru...	AA								
ABCDEF	Image Instru...	ABCDEF								
R14271/07	Image Instru...	ABELS								
18630	Image Instru...	ANONYMO								
3M_JBT	Image Instru...	APC								
21150	Image Instru...	ARAT								
10103	Image Instru...	ATRON								
Aachen	Image Instru...	Aachen								
270682	Image Instru...	Aalborg								
0705	Image Instru...	Techno Techno								

From:

<https://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

https://www.onyxwiki.net/doku.php?id=program_onyxstat&rev=1587813728

Last update: 2020/04/25 13:22

