




# WebViewer Standard Controls and Examples




Webview Example Aligner planning






















Webview Example Facescan



Webview Example Headscan



Icon	Selection / Alignment	Icon	Animation
	Display upper object group / Bottom view		Start animation
	Display entire object / Frontal view		Stop animation
	Display lower object group / Top view		Play animation until next step
	Display rotation Frontal-Right-Rear-Left-Bottom-Top	<b>Icon</b>	<b>Speed</b>
<b>Icon</b>	<b>View</b>		Play animation /slow
	Display vertex colors: on / offs		Play animation /medium
	Display attachments: off / transparent / on		Play animation /fast
	Display malocclusion: off / transparent / on	<b>Icon</b>	<b>Loop</b>
	Display: step name / step name + comment / off		Repeat animation when finished
	Display ruler: on / off		Stop animation when finished
	Clipping along Y: none / right / left	<b>Icon</b>	<b>Info</b>
			Display available patient info for the displayed case

#### Time line



Step 1-3 upper only - Step 4-6 upper/lower - Step 7-8 lower only - Step 9-13 upper/lower

## Links



From:  
<http://onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:  
[http://onyxwiki.net/doku.php?id=en:3d\\_viewn](http://onyxwiki.net/doku.php?id=en:3d_viewn)

Last update: **2023/03/02 15:46**

