


## #4314 Scalable Base Trays

Since release 3.2.78, in module [Cast Adjust 3D](#), all selectable standard base tray types can additionally be downscaled:

- 100%
- 90%
- 80%



This should allow a better adjustment between base tray and individual dental arch size and therefore to save material for rapid prototyping of study models.

From:

<https://onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

<https://onyxwiki.net/doku.php?id=en:4314>

Last update: 2017/06/14 14:10

