


Module Add Image 3D

In Module Add Image 3D different types of 3D data sets can be imported from pre-defined image sources, classified acc. to [session](#) and [image type](#) and completed with [attributes and comments](#) before saved to the Onyx database.

First Steps

1. Activate "Clipping (Line)" in tab "Clip Object"
2. Place Landmarks at the model border, last point close to the first one
3. Press button "Trim" and confirm the correct trimming
4. Select and add socket in tab "Model Base"
5. Save as new finding

Module user interface



Links



https://image-instruments.de/Onyx3/Handout_Digitize_MA_EN.pdf



<https://youtu.be/TAAGkmFIKD4>



<https://youtu.be/Q9R9V5QgEoE>



[Import | Adjust Scan](#)

From:

<http://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

http://www.onyxwiki.net/doku.php?id=en:add_image_3d&rev=1585253459

Last update: **2020/03/26 21:10**

