


Module Combine 3D




Program version	BASIC	2D PRO	3D PRO	LAB	OMS
Availability					

 not included

 included by

default

 optionally


included

The Combine 3D module provides 2 separate modes for the combination of two 3D findings based on overlapping features to allow additional options for visualization, diagnostics or planning.



First Steps

1. Open module with 2 findings
2. Select corresponding reference points
3. For ICP registration, select surface region
4. Switch to tab |Combine| and configure result
5. Save result as finding

Module user interface



Verweise

	Combine 3D
	Panel Surface Registration 3D

From:
<http://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
http://www.onyxwiki.net/doku.php?id=en:combine_3d&rev=1667208413

Last update: **2022/10/31 10:26**

