

Activation

◀ TP OnyxAdmin



Note: This instruction is only needed if a new clinic server has to be installed

1. Import Activation Key

To activate the license with the activation key received by email, copy the key to clipboard (CTRL+C) and paste by button [Clipboard] in Onyxregister step 3.

	License request Step 3
	License activation - Animation

Note that the initial activation can only be completed if application OnyxCeph³(TM) was locally closed before. OnyxRegister is a separate executable and will stay running.

From:
<https://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
https://onyxwiki.net/doku.php?id=en:moh:tp_admin4

Last update: **2022/03/10 08:40**

