

Client Application OnyxRegister



Client application OnyxRegister is used to [activate or renew a license](#). The program is opened via [main menu](#) [Options | Register or Options | Renew license], [info menu](#) [License] or from program group Onyx3Ceph in the Windows™ start menu at the local workstation.

User interface

Image Instruments

OnyxCeph³™ Registration

Step 1 Step 2 Step 3

Save and Send Registration File

Checksum: 4280

Practice Name: Blue Tooth Ortho Clinic

Installation Code	Version	Server Name	Release
61JL-M7J7	73	DB1	3.2

Title: Dr. Phone 1: +666-951753

First Name: Blue Phone 2:

Last Name: Tooth Fax: 0371 9093149

Address: 123 Main street Email: info@r666.digital

Zip Code: PA, 17101 Internet: www.r666.digital

City: Anycity VAT - ID:

Country: UK Payment: Bank Transfer

Version: Lab License: OSL

Modules: VTO3D, Waefer, BondTray, WireBonding, Aligner, Retainer, Sim3D, TADmatch, BracketErase, OrthoApps, BiteSplint

For registration save the information below to file or to clipboard. Send the registration file or the content of the clipboard by email to:

register@onyx-ceph.de

Email... Save... Copy...

Help Close

Links

 [License Activation](#)

From:
<http://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
http://www.onyxwiki.net/doku.php?id=en:program_onyxregister

Last update: **2023/03/02 17:22**

