

Module Retainer 3D



Program version	BASIC	2D PRO	3D PRO	LAB	OMS
Availability	✗	✗	📦	📦	📦

- ✗

enthalten

nicht
- 📦

enthalten

immer
- 📦

enthalten


optional

Module Retainer 3D module can be used to construct virtual retainers based on an unsegmented or setup model created with the V.T.O.3D module.




First Steps

1. Place three points defining the front point and left/right end of the retainer
2. Set parameters such as height, width and precision
3. Compute the retainer for each jaw
4. Export the retainer as 2D path, 3D object or send it to a service provider
5. Save project

Module user interface



Links

	Der krönende Abschluss der KFO-Behandlung
	Module Retainer
	YOAT Benderl-Interface

From:
<http://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
http://www.onyxwiki.net/doku.php?id=en:retainer_3d&rev=1644236892

Last update: **2022/02/07 13:28**

