

# Module Retainer 3D



Program version	BASIC	2D PRO	3D PRO	LAB	OMS
Availability	✗	✗	📦	📦	📦

✗

not included

📦

included by

📦

optionally

📦


included

Module Retainer 3D can be used to construct virtual retainers based on an unsegmented or setup model created with the V.T.O.3D module.

## First Steps

1. Place three points defining the front point and left/right end of the retainer
2. Set parameters such as height, width and precision
3. Compute the retainer for each jaw
4. Export the retainer as 2D path, 3D object or send it to a service provider
5. Save project

## Module user interface



## Links

	<a href="#">Der krönende Abschluss der KFO-Behandlung (in German)</a>
	<a href="#">Module Retainer</a>
	<a href="#">YOAT Bender-I-Interface</a>
	<a href="#">Digital Retainer Planung: Software Requirements</a>

From:

<http://www.onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

[http://www.onyxwiki.net/doku.php?id=en:retainer\\_3d&rev=1676910771](http://www.onyxwiki.net/doku.php?id=en:retainer_3d&rev=1676910771)

Last update: **2023/02/20 17:32**

