

Use OnyxReport

Patient

FB2354 - Roots, Tooth
4178 - Klüber, Johannes

Daten

☒ Send Anonymously

Options

Recipient's language
English

☐ Password

☒ Send Email with FTP-Upload

Images **Slides**

No.	Type	Date of Record	Image
<input checked="" type="checkbox"/> 1	CT-HeadScan	03.03.2022	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/> 2	Lateral Right	03.11.2021	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/> 3	Cast Deciduous 3D	23.11.2021	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/> 4	Cast Mixed 3D	23.11.2021	<input checked="" type="checkbox"/>

Comment about Patient

Add patient specific comments here

Sender

Image Instruments

Message for recipient

Add a message for the recipient here


SAMPERMANS IOS INTERNATIONAL ORTHODONTIC SEMINARS

Help Preview Close

Controls

- [Auswahlpanel](#)
- [Iconmenü](#)
- [Aufruf einer Animation](#)

Selection Panel




Icon Menu




Calling an Animation

Please start the full screen of the selected finding first. By clicking the right mouse button you can now reach the context menu 'Animation' (see below).



In the panel, the changes can be displayed using the scroll bar or the animation can be started using the button.



Links

	Client Application OnyxReport
	OnyxReport

From:

<http://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

http://www.onyxwiki.net/doku.php?id=en:use_onyxreport

Last update: **2022/04/27 08:00**

