

# Systemoptionen - Tab |Daten hochladen|

Auf Systemoptionen-Reiter |Daten hochladen| erfolgen alle Einstellungen für die Kommunikation mit Webservern über FTP- und HTTP-Protokoll zum Hoch- bzw. Herunterladen von [Container](#)-, Report- und [Webviewer](#)-Dateien.

The screenshot shows the 'System Options' dialog box with the 'Container' tab selected. The left sidebar contains a list of categories: Patient Data, Image Types, Sessions, Galleries, Printouts, Forms, Display, Interfaces, Third Party Tools, Event, Upload Data (highlighted), and Email. The main area of the dialog has four tabs: Container, Report, Import, and 3D-Web-Viewer. The 'Container' tab is active, showing fields for Host, Connection Security (set to None), User Name, User Password (masked with dots), Host Directory, and Passive Mode (checked). Below these fields is a section for 'Download - Protocol / Data Address' with a dropdown menu set to 'https://' and a text field containing 'onyx.orthorobot.com/container'. At the bottom of the dialog are 'Help' and 'Close' buttons.

Category	Field	Value
Container	Host:	
	Connection Security:	None
	User Name:	
	User Password:	.....
	Host Directory:	
	Passive Mode:	<input checked="" type="checkbox"/>
Download - Protocol / Data Address:	Protocol	https://
	Data Address	onyx.orthorobot.com/container

[Click auf Reiter öffnet zugehörigen Wiki-Seite]

From:

<http://www.onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

[http://www.onyxwiki.net/doku.php?id=tab\\_uploaddata&rev=1617877397](http://www.onyxwiki.net/doku.php?id=tab_uploaddata&rev=1617877397)

Last update: **2021/04/08 12:23**

