


## FAQ: Display 3D Dataset in Original Size

### Question:

How can I display a 3D object in original size on the monitor?

### Answer:

Right click onto the mid position of the right zoom slider in each 3D view panel:



Monitor calibration can be adjusted by menu Options|System Optionen on Tab |Digitize|.

From:

<https://www.onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

<https://www.onyxwiki.net/doku.php?id=en:faq10010>

Last update: **2018/03/16 11:27**

