OnyxDBServer

The OnyxDBServer program canbe used to make settings for database management and data backup. It is set up with the server installation and, if possible, should only be started or stopped as a service or program using the Onyx3Config program.

It is recommended to run the program as servvice to ensure that it can be used independently of an active Windows[™] user login. However, settings can only be made or changed if OnyxDBServer was started as a program. In such case, the main application window is minimized in the system tray of the task bar and can be started from there via the context menu.

Program user interface

OnyxCeph ³ ™ - Database Server			-		
rver User Databases Status	s Settings	Messages	Backup		
itatus Information					
OnyxCeph³™ -	Database S	Server (ac	tive)		
Server Version:	3.2.12 (45)			
Database Version:	1				
Server Name:	[RK8]				
Server Time:	05.04.20	20 14:30:40	D		
Server Running Time:	0.01:04:	50			
TCP/IP - Network Address:	172.16.1	7.25			
'Named Pipe' - Network Name:	RK8				
Active Sessions:	0				
Stop Server	C Refresh		0	Close	e

Server settings can be made or adjusted on the following tab sheets:

- |Server|
- |User|
- |database|
- |Status|
- |Settings|
- |Messages|
- |Backup|

From: https://www.onyxwiki.net/ - [OnyxCeph³™ Wiki]

Permanent link: https://www.onyxwiki.net/doku.php?id=en:program_onyxdbserver&rev=1586091603

Last update: 2020/04/05 15:00

