


Client Application OnyxView



Client application OnyxView provides a thumbnail overview of all 2D and 3D image findings available for the selected patient in a session view on separate tabs. A single image viewer can be opened in a separate window for the thumbnails selected there.

OnyxView can be opened from all module windows using the icon button. The single image viewer can also be used in the main window on the tab | Images / 2D data | and | Images / 3D data | via thumbnail context menu [Send to | Image Viewer].

Thumbnail-Fenster



From:
<http://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
http://www.onyxwiki.net/doku.php?id=en:program_onyxview

Last update: 2023/02/16 22:24



