



This page is not fully translated, yet. Please help completing the translation.
 (remove this paragraph once the translation is finished)


Module Ricketts V.T.O.

Module Ricketts V.T.O. [according to the systematics developed by R.M. RICKETTS] can be used to to simulate growth and certain orthodontic and/or surgical treatment measures with the help of a lateral ceph tracing in combination with an associated lateral photo in their effect on cephalometric measurements, space discrepancy and soft tissue profile.

First Steps

1. Modul mit FRS und (optional) zugehörigem Profilbild öffnen
2. Vorhersagezeitraum und Wachstumsrichtung festlegen
3. Behandlungsmaßnahmen simulieren
4. Änderungen der Durchzeichnung analysieren (grafisch, numerisch)
5. Weichteilsimulation starten: Vorher-Nachher-Vergleich
6. Planung als Projekt speichern
7. Modifizierte Befunde (FRS, Foto) speichern

Module user interface



Links



From:

<https://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

https://www.onyxwiki.net/doku.php?id=en:ricketts_vto_2d&rev=1584543144

Last update: **2020/03/18 15:52**

