


Interface Configuration Digital Xray: VixWin Scan



Parameter	Setting
[Support Activated]	Activates import icon button in module Add Image 2D
Programmpfad [...]	Lokaler Pfad zur ausführbaren Datei ..\VixWin.exe
[Kopi speichern]	Bildkopie in VixWin speichern
[Bildformat]	Auswahl Bildspeicherformat

From:
<http://onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
http://onyxwiki.net/doku.php?id=en:vixwin_scan

Last update: **2022/06/27 07:38**

