

Client-Programm OnyxRegister



Das Client-Programm OnyxRegister dient der [|Aktivierung der Lizenz](#). Der Programmaufruf erfolgt über das [Hauptmenü](#) [Optionen|Registrieren bzw. Optionen|Lizenz erneuern], das [Infomenü](#) [Lizenz] oder aus der Programmgruppe OnyxCeph^{3™} im Windows-Startmenü am lokalen Arbeitsplatz.

Bedienoberfläche

Image Instruments

OnyxCeph³™ Registration

Step 1 Step 2 Step 3

Save and Send Registration File

Checksum: 4280

Practice Name: Blue Tooth Ortho Clinic

Installation Code: 61JL-M7J7 Version: 73 Server Name: DB1 Release: 3.2

Title: Dr. Phone 1: +666-951753

First Name: Blue Phone 2:

Last Name: Tooth Fax: 0371 9093149

Address: 123 Main street Email: info@r666.digital

Zip Code: PA, 17101 Internet: www.r666.digital

City: Anycity VAT - ID:

Country: UK Payment: Bank Transfer

Version: Lab License: OSL

Modules: VTO3D, Waefer, BondTray, WireBonding, Aligner, Retainer, Sim3D, TADmatch, BracketErase, OrthoApps, BiteSplint

For registration save the information below to file or to clipboard. Send the registration file or the content of the clipboard by email to:

register@onyx-ceph.de

Email... Save... Copy...

Help Close

Verweise

[Lizenz-Aktivierung](#)

From:
<http://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
http://www.onyxwiki.net/doku.php?id=program_onyxregister

Last update: **2023/03/02 16:32**

