

# OnyxStat

Mit Hilfe des Dienstprogramms OnyxStat können Onyx-Datenbanken statistisch nach unterschiedlichen Kriterien durchsucht und ausgewertet werden. Die Abfrageergebnisse können gespeichert und im CSV-Format exportiert werden.

OnyxStat erfordert eine separate, einmalig kostenpflichtige Feischaltung.

Der Programmstart erfolgt aus der Programmgruppe OnyxCeph<sup>3</sup>™ im Windows-Startmenü am lokalen Arbeitsplatz.

The screenshot displays the OnyxStat application window. It features a menu bar (Project, Evaluation, Query, Window, Options, Help) and a toolbar. The main area is divided into two panes. The left pane, titled 'Query 1', shows a table with columns: Patient ID, Patient Client, Patient Name, Patient First Name, Patient Gender, Patient Birthday, Patient Attributes, Finding Date, Image Type, 2D/3D, Attributes, and Status. The right pane, titled 'Query 2', shows a similar table with columns: Patient ID, Client, Name, First Name, Gender, Birth, Finding Date, Image Type, 2D/3D, Attributes, and Status. A third, smaller window is overlaid on the 'Query 2' table, showing a detailed view of a specific patient's data. On the right side of the application, there is a 'View' panel with options for 'Color', 'Categories as Text', 'Sort By', 'Patient Filter', 'Data Filter', and 'Columns'. Below these options is a 'Go to' section with a dropdown menu set to 'Image', and a small thumbnail image of a patient's head scan.

From:

<https://www.onyxwiki.net/> - [OnyxCeph<sup>3</sup>™ Wiki]

Permanent link:

[https://www.onyxwiki.net/doku.php?id=program\\_onyxstat&rev=1515349261](https://www.onyxwiki.net/doku.php?id=program_onyxstat&rev=1515349261)

Last update: **2018/01/07 19:21**

