

OnyxStat



Mit Hilfe des Dienstprogramms OnyxStat können Onyx-Datenbanken statistisch nach unterschiedlichen Kriterien durchsucht und ausgewertet werden. Die Abfrageergebnisse können gespeichert und im CSV-Format exportiert werden.

OnyxStat erfordert eine separate, einmalig kostenpflichtige Feischaltung.

Der Programmstart erfolgt aus der Programmgruppe OnyxCeph3 im Windows-Startmenü am lokalen Arbeitsplatz.

Bedienoberfläche

The screenshot displays the OnyxStat application window. The main area shows a table with patient data. A 'Query 2' window is open, displaying a filtered view of the data. The sidebar on the right contains a 'Query 2' section with a 'View' button, a 'Color' button, and a 'Categories as Text' checkbox. Below this are 'Sort By', 'Patient Filter', 'Data Filter', and 'Columns' sections. The 'Columns' section lists various fields with checkboxes, including Patient ID, Client, Name, First Name, Gender, Birth, Ceph, 2D, Attributes, Status, Finding Date, Image Type, 2D/3D, and Age. At the bottom of the sidebar is a 'Go to' section and an 'Image' section showing a patient's lateral cephalogram.

Patient ID	Client	Name	First Name	Gender	Birth	Ceph	2D	Attributes	Status
70020	Image Instru...	1	2	Male	24.04.2011				Initial
7654321	Image Instru...	1		Male	20.03.2003				
10015	Image Instru...	15							
10017	Image Instru...	17							
54321	Image Instru...	2							
10091	Image Instru...	2							
2000	Image Instru...	2D Ceph A							
ID001344	Image Instru...	3.7.1117.							
70016	Image Instru...	3D							
3Q	Image Instru...	3Q Techno							
10128	Image Instru...	3d							
10178	Image Instru...	3shape							
10118	Image Instru...	<SkX							
10115	Image Instru...	A							
AAAAA	Image Instru...	AA							
AAAAA	Image Instru...	AA							
AAAAA	Image Instru...	AA							
ABCEDEFG	Image Instru...	ABCEDEFG							
R14271/07	Image Instru...	ABELS							
18630	Image Instru...	ANONYMO							
3M_IBT	Image Instru...	APC							
21150	Image Instru...	ARAT							
10103	Image Instru...	ATRON							
Aachen	Image Instru...	Aachen							
270682	Image Instru...	Aalborg							
0706	Image Instru...	Aalborg							

From:

<https://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:

https://www.onyxwiki.net/doku.php?id=program_onyxstat&rev=1587923555

Last update: **2020/04/26 19:52**

