

OnyxStat



Mit Hilfe des Dienstprogramms OnyxStat können Onyx-Datenbanken statistisch nach unterschiedlichen Kriterien durchsucht und ausgewertet werden. Die Abfrageergebnisse können gespeichert und im CSV-Format exportiert werden.

OnyxStat erfordert eine separate, einmalig kostenpflichtige Feischaltung.

Der Programmstart erfolgt aus der Programmgruppe OnyxCeph3 im Windows-Startmenü am lokalen Arbeitsplatz.

Bedienoberfläche

The screenshot shows the OnyxStat application window. The main data table has columns: Patient ID, Patient Client, Patient Name, Patient First Name, Patient Gender, Patient Birthday, Patient Attributes, Finding Date, Image Type, 2D/3D, Attributes, and Status. A 'Query 2' dialog box is overlaid, showing a filtered table with columns: Patient ID, Client, Name, First Name, Gender, Birth, Finding Date, Image Type, 2D/3D, Attributes, and Status. The 'View' panel on the right includes a 'Color' button, a 'Categories as Text' checkbox, a 'Sort By' dropdown, and a 'Patient Filter' dropdown. The 'Columns' section lists various fields with checkboxes. The 'Go to' section has a radio button for 'Image' and a small image thumbnail.

From:
<https://www.onyxwiki.net/> - [OnyxCeph³™ Wiki]

Permanent link:
https://www.onyxwiki.net/doku.php?id=program_onyxstat&rev=1588053732

Last update: 2020/04/28 08:02



